

codelink

November 11, 2009

R topics documented:

arrayNew	2
arraySize	2
as.matrix.Codelink	3
bkgdCorrect	3
Codelink-class	4
codelink.example	5
codelink.exprset	6
createWeights	6
cutCV	7
decDetect	7
dim.Codelink	8
fc2Cond	9
imageCodelink	9
logCodelink	10
mergeArray	11
na2false	12
Nondocumented-objects	12
normalize.loess	13
normalize	14
plotCorrelation	15
plotCV	15
plotDensities	16
plotMA	17
printHead	18
readCodelink	19
readHeader	20
readHeaderXLS	20
reportCodelink	21
selCV	22
SNR	22
writeCodelink	23

Index	24
--------------	-----------

`arrayNew`*Create a new x11 device*

Description

Create a new x11 device with dimensions suited to be used with `imageCodelink()`.

Usage

```
arrayNew(f=2, chip="rwgcod")
```

Arguments

<code>chip</code>	character; Codelink chip to be used in <code>imageCodelink</code> .
<code>f</code>	numerical; scaling factor.

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  arrayNew()  
  imageCodelink(codelink.example)  
## End (Not run)
```

`arraySize`*Determine the size of the array*

Description

When loading Codelink arrays in text format (as exported from the Codelink software) this function retrieves the correct size of the array. This is useful because those files contain an indetermined number of empty lines at the end. Thus, reading the entire data matrix doesn't work.

Note

Not meant to be used directly.

Author(s)

Diego Diez

as.matrix.Codelink *Return a matrix of intensity values*

Description

Takes a Codelink object and returns a matrix with the intensity values available.

Usage

```
## S3 method for class 'Codelink':  
as.matrix(x, ...)
```

Arguments

x an object of class "Codelink".
... additional arguments added to generic as.matrix since R-2.5.x

Value

A matrix with the intensity values.

Author(s)

Diego Diez

See Also

[as.matrix](#)

Examples

```
## Not run:  
  data(codelink.example)  
  mat <- as.matrix(codelink.example)  
  is(mat)  
## End(Not run)
```

bkgdCorrect *Background correction of intensity values.*

Description

Takes a Codelink object with Spot mean and Bkgd median values and corrects applying one of the methods available.

Usage

```
bkgdCorrect(object, method = "half", preserve = FALSE, verbose = FALSE,  
          offset = 0)
```

Arguments

<code>object</code>	character; an object of class "Codelink".
<code>method</code>	character; the correction method to use, one of "none", "subtract", "half" and "normexp".
<code>preserve</code>	logical; if Smean and Bmedian slots should be preserved.
<code>verbose</code>	logical; if TRUE print some information with method normexp.
<code>offset</code>	numeric; value to add to intensities.

Details

Available methods are: . none: left intensities untouched. . subtract: simple subtraction of Bkgd median from Spot mean. . half: the same as above but avoid negative values setting all intensity values below zero to 0.5. . normexp: apply normexp background adjustment from package limma.

Value

An object of class Codelink with corrected intensity values, that is Ri slot.

Author(s)

Diego Diez

Examples

```
## Not run:
  data(codelink.example)
  codelink.example <- bkgdCorrect(codelink.example, method = "half")
## End(Not run)
```

Codelink-class *Class Codelink*

Description

This is the storage class for Codelink data

Objects from the Class

Object are created after reading text codelink files with read.Codelink()

Slots

sample: Object of class "character" containing the sample names
file: Object of class "character" containing the file names
name: Object of class "character" containing the probe ID
method: Object of class "list" containing log information
Smean: Object of class "matrix" containing spot mean intensities
Bmedian: Object of class "matrix" containing background median intensities
Ri: Object of class "matrix" containing raw intensities

ni: Object of class "matrix" containing normalized intensities

snr: Object of class "matrix" containing signal to noise ratio values

flag: Object of class "character" containing assigned flags

Methods

Note

More details are contained in the vignette

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
## End(Not run)
```

codelink.example *Dataset of class 'Codelink'*

Description

Dataset from a h20kcod (Codelink Human UniSet I 20k) array containing 2 samples and ~20000 probes.

Usage

```
data(codelink.example)
```

Format

A [Codelink](#) object containing 20469 probes and 2 samples.

Author(s)

Diego Diez

`codelink.exprset` *Dataset of class 'CodelinkSet'*

Description

Dataset from a h20kcod (Codelink Human UniSet I 20k) array containing 4 samples and ~20000 probes.

Usage

```
data(codelink.example)
```

Format

A `Codelink` object containing 20469 probes and 4 samples.

Author(s)

Diego Diez

`createWeights` *Create weight for Codelink chips*

Description

Create a weight matrix to be used in limma.

Usage

```
createWeights(object, type=NULL)
```

Arguments

<code>object</code>	an object of class "Codelink".
<code>type</code>	weight assigned to each Probe_type.

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  w <- createWeights(codelink.example, type = list(FIDUCIAL = 0.1))  
## End(Not run)
```

cutCV	<i>Calculate cutoff based in CV.</i>
-------	--------------------------------------

Description

Takes a Codelink object and calculate cutoff based in CV.

Usage

```
cutCV(object, subset=c(1:dim(object)[2]))
```

Arguments

object	an object of class "Codelink".
subset	subset of arrays to calculate cutoff with

Details

First it computes the median of CV for each gene over all arrays. Then it computes the mean and sd of all medians. Finally:

```
cutoff = mean + 3 * sd
```

Author(s)

Diego Diez

Examples

```
## Not run:  
# data: Normalized Codelink object merged.  
cutoff <- cutCV(data)  
## End(Not run)
```

decDetect	<i>Determine decimal type of Codelink files</i>
-----------	---

Description

Determine decimal type of Codelink files.

Usage

```
decDetect(file, nlines)
```

Arguments

file	the file to be read.
nlines	number of lines to skip.

Value

Decimal type.

Author(s)

Diego Diez

dim.Codelink

Return the dimension of a Codelink object.

Description

Takes a Codelink object and returns the dimension (genes x samples).

Usage

```
## S3 method for class 'Codelink':  
dim(x)
```

Arguments

x an object of class "Codelink".

Value

A numeric vector with the dimensions.

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  dim(codelink.example)  
## End (Not run)
```

fc2Cond	<i>Select probes based on fold change calculation</i>
---------	---

Description

Takes a Codelink object and calculate fold changes (M) between two conditions (samples). Then select genes based on those who pass the passed cutoff.

Usage

```
fc2Cond(object, cond1=NULL, cond2=NULL, fc=1.0, verbose=FALSE)
```

Arguments

object	an object of class "Codelink".
cond1	numeric or character; First condition to compute M.
cond2	numeric or character; Second condition to compute M.
fc	value of the fold change cutoff
verbose	logical; if some information is printed on the console.

Details

Conditions can be passed as characters or as numeric index from the sample slot. The intensities are internally transformed to log2 if needed. The M value is computed as:

$$M = \text{cond1} - \text{cond2}$$

Value

A logical vector indicating which genes pass the cutoff

Author(s)

Diego Diez

imageCodelink	<i>Image plot of Codelink arrays</i>
---------------	--------------------------------------

Description

Plot and image of a Codelink array if the layout information is found.

Usage

```
imageCodelink(object, array = 1, what = "bg",
  low="black", high="white", mar=c(1,1,1,1),
  gr=1, gc=1, log.it=FALSE, ...)
```

Arguments

object	an object of class "Codelink".
array	array to be used.
what	with data plot: bg, smean, ri, ni.
low	color used for low intensities.
high	color used for high intensities-
mar	character vector specifying margins.
gc	numerical; number of grid columns.
gr	numerical; number of grid rows.
log.it	logical; if TRUE data is log2 transformed (if not yet).
...	additional arguments passed to image.

Author(s)

Diego Diez

Examples

```
## Not run:
  data(codelink.example)
  imageCodelink(codelink.example)
## End(Not run)
```

logCodelink

Application of logCodelink to Codelink object

Description

Takes a Codelink object and apply logCodelink to intensity values.

Usage

```
logCodelink(object)
```

Arguments

object	an object of class "Codelink" or a list of genes.
--------	---

Value

A Codelink object with logCodelink intensities.

Author(s)

Diego Diez

See Also

[log2](#)

Examples

```
## Not run:
  data(codelink.example)
  codelink.example <- logCodelink(codelink.example)
## End(Not run)
```

mergeArray	<i>Merge Codelink Bioarrays Data</i>
------------	--------------------------------------

Description

Merge data in a Codelink Object corresponding to same samples. Need a vector indicating the classes and an optional vector indicating the labels of the merged samples.

Usage

```
mergeArray(object, class, names=NULL, method="mean",
           log.it=FALSE, merge.snr=TRUE)
```

Arguments

object	an object of class "Codelink".
class	a numeric vector indicating the classes.
names	an optional character vector indicating labels for each class.
method	the method used to summarize. Currently only "mean" supported.
log.it	logical; a logical indicating if log ₂ values should be returned.
merge.snr	logical; a logical indicating if SNR values should be merged.

Value

An object of class "Codelink".

Author(s)

Diego Diez

Examples

```
## Not run:
  data(codelink.example)
  codelink.example <- bkgdCorrect(codelink.example)
  codelink.example <- normalize(codelink.example, log.it = FALSE)
  codelink.example <- mergeArray(codelink.example, class = c(1,1),
                                names = "SAMPLE", log.it = TRUE)
## End(Not run)
```

`na2false`*Set NAs to FALSE*

Description

Takes a logical vector as input and set all NAs to FALSE. This may happens when comparison is done on NA values.

Usage

```
na2false(x)
```

Arguments

`x` a logical vector.

Value

A logical vector without NAs.

Author(s)

Diego Diez

Examples

```
## Not run:
a <- c(1, 2, 3, NA, 5)
b <- c(5, 4, NA, 2, 1)
sel <- a > b
sel <- na2false(sel)
## End(Not run)
```

`Nondocumented-objects`*Here goes every undocumented function.*

Description

Every function without help page is redirected here. Either because it is not intended to have a help page or because it is a bug. If you think it is the later please contact the package maintainer.

NOTE: Actually it is none of that reasons but the fact that the package is under heavy development.

Author(s)

Diego Diez

normalize.loess *Normalized chips using cyclic loess.*

Description

Takes a matrix and apply cyclic loess normalization. It is based in normalize.loess from package affy but supports NA.

Usage

```
normalize.loess(mat, subset = sample(1:(dim(mat)[1]),
  min(c(5000, nrow(mat)))), epsilon = 10^-2, maxit = 1, log.it = TRUE,
  verbose = FALSE, span = 2/3, family.loess = "symmetric", weights = NULL)
```

Arguments

mat	a matrix with columns containing the values of the chips to normalize.
subset	a subset of the data to fit a loess to.
epsilon	a tolerance value (supposed to be a small value - used as a stopping criterium).
maxit	maximum number of iterations.
log.it	logical. If TRUE it takes the log2 of mat
verbose	logical. If TRUE displays current pair of chip being worked on.
span	parameter to be passed the function <code>loess</code>
family.loess	parameter to be passed the function <code>loess</code> . "gaussian" or "symmetric" are acceptable values for this parameter.
weights	a vector of weights for the individual measurements.

Value

A matrix of normalized values.

Author(s)

Diego Diez

Examples

```
## Not run:
  mat <- matrix(sample(500), 100, 5)
  mat <- normalize.loess(mat)
## End(Not run)
```

`normalize`*Normalization wrapper for Codelink objects.*

Description

Takes a Codelink object and applies normalization to intensity values.

Usage

```
normalize(object, method = "quantiles", log.it = TRUE, preserve = FALSE,
         weights = NULL, verbose = FALSE)
```

Arguments

<code>object</code>	an object of class "Codelink".
<code>method</code>	method to use in normalization.
<code>log.it</code>	logical; if data should be log2.
<code>preserve</code>	logical; if Ri slot should be preserved.
<code>weights</code>	weights vector for method CyclicLoess.
<code>verbose</code>	should informative output be printed.

Details

Currently supported methods include "loess", "quantiles" and "median". Median normalization is analogous to the default method applied for the manufacturer in the Codelink software. Loess is a modified version of CyclicLoess implemented in the affy package, allowing missing values and weights. Quantile normalization uses the `normalizeQuantiles()` function in the limma package.

Value

A Codelink object with normalized intensity values.

Author(s)

Diego Diez

Examples

```
## Not run:
  data(codelink.example)
  # Background correction.
  codelink.example <- bkgdCorrect(codelink.example,
  method = "half")
  # Normalization.
  codelink.example <- normalize(codelink.example,
  method = "quantile")
## End(Not run)
```

plotCorrelation *Plot correlation scatterplot between two arrays*

Description

Takes a Codelink object as argument and plot Correlation scatterplot of two arrays.

Usage

```
plotCorrelation(object, x=1, y=2, cutoff=FALSE, label="type", title=NULL)
```

Arguments

object	an object of class "Codelink".
x	array to be used in x axis.
y	array to be used in y axis.
cutoff	cutoff used to show fold change.
label	labels to shown.
title	The title of the plot.
xlim	range for the X axis.
ylim	range for the Y axis.

Author(s)

Diego Diez

See Also

[plot](#)

Examples

```
## Not run:
  data(codelink.example)
  plotCorrelation(codelink.example)
## End(Not run)
```

plotCV *Plot of CV*

Description

Takes a Codelink object and plot de distribution of CV after applying mergeCodelink.

Usage

```
plotCV(object, subset=c(1:dim(object)[2]), cutoff=NULL, title=NULL, lege
```

Arguments

object	an object of class "Codelink".
subset	subset of arrays to plot
cutoff	cutoff of CV to be shown.
title	title of the plot.
legend.cex	factor to apply to the fonts in the legend to fit.

Author(s)

Diego Diez

plotDensities *Plot Densities*

Description

Takes a Codelink object and plot the distributions of intensities.

Usage

```
plotDensities(object, subset=1:dim(object)[2], title=NULL,  
              legend.cex=1, what=NULL)
```

Arguments

object	an object of class "Codelink".
subset	subset of arrays to be plotted (default: all).
title	title of the plot.
legend.cex	font factor use in legend to fit.
what	what data to plot, may be "bg", "smean", "snr", "ri" or "ni"

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  plotDensities(codelink.example)  
## End(Not run)
```


plotMA

*MA plot***Description**

Takes a Codelink object and plot M vs A.

Usage

```
plotMA(object, array1 = 1, array2 = NULL, cutoff = c(-1, 1), label = NULL,
       type = NULL, high.list = NULL, high.col = "blue", high.pch = 21,
       high.bg = "cyan", snr = NULL, snr.cutoff = 1, legend.x = NULL, pch
       subset = NULL, title = NULL, xlim = NULL, ylim = NULL)
```

Arguments

object	an object of class "Codelink" or "MAarrayLM".
array1	first array to be used.
array2	second array to be used.
cutoff	cutoff to be used as fold change marker.
label	type of labeling used in legend.
type	spot type information.
high.list	list of genes highlighted.
high.col	color used for high genes.
high.pch	pch used for high genes.
high.bg	background color used for high genes.
snr	vector with SNR values, usually, taking rowMeans() from a SNR matrix.
snr.cutoff	SNR cutoff used for label spots.
legend.x	relative position of the legend.
pch	pch style used to main spots.
subset	subset of spots used to plot based on 'type' slot.
title	title of the plot.
xlim	range for the X axis.
ylim	range for the Y axis.

Details

This function has suffered recent re-working, to increase the usability and to clean a little bit the code.

If array2 is NULL a median array is computed using all available arrays. Then the values of M and A are computed using the following formula:

$$M = \text{array2} - \text{array1}$$

$$A = (\text{array2} + \text{array1}) / 2$$

If type information is available in the Codelink object, or provided through the 'type' argument, spots are colored based on that. DISCOVERY spots are plotted black with pch = "." whereas the

other classes are plotted with different background colors, using gray as border to increase contrasts. For that `pch = 21` is used. If `snr` is specified as label option, the SNR is used to label spots, if available in the Codelink object. In this case, the mean SNR across all arrays is used when `array2 = NULL`.

Some parameters may not be working right now, as the new function is using a different method to labels spots.

The legend is 'automagically' located, but this can be overridden with the `legend.x` argument.

In addition, a subset of the spots can be plotted based on type information when available. This allows, for example, to plot only DISCOVERY spots.

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  plotMA(codelink.example)  
## End(Not run)
```

`printHead`

Print briefly a Codelink object

Description

Takes a Codelink object and print a summary information of the data estored. It is based on `printHead()` from package `limma`.

Usage

```
printHead(x)
```

Arguments

`x` an object of class "Codelink".

Author(s)

Diego Diez

readCodelink	<i>Read Codelink Bioarrays Data</i>
--------------	-------------------------------------

Description

Read data exported as text by Codelink Software. It reads values (normalized by Codelink Software or not) flags and information about probes.

Usage

```
readCodelink(files=list.files(pattern="TXT"), sample.name=NULL,
             flag, dec=NULL, type="Spot", preserve=FALSE, verbose=2,
             file.type="Codelink", check=TRUE, fix=FALSE)
```

Arguments

files	list of files to read.
sample.name	vector of same length as files with sample names.
flag	list with values to assign based on Flag quality values.
dec	character indicating the decimal character used in the files.
type	character indicating which base value to read from files.
preserve	logical, if TRUE Bkgd_stdev slot is not removed (if present).
verbose	numerical, set the level of information. Level 2 set as old behaviour. Level > 2 output some debug info.
file.type	exported file type, currently Codelink or XLS file formats supported.
check	logical, check for probe order consistency.
fix	logical, try to fix probe order consistency.

Value

An object of class "Codelink".

Author(s)

Diego Diez

See Also

[read.table](#)

Examples

```
## Not run:
# reading default extension (TXT).
data <- readCodelink()

# specify a different one.
files <- list.files(pattern = "txt")
data <- readCodelink(files = files)
```

```
# example.
data(codelink.example)
## End(Not run)
```

readHeader *Read Header from Codelink Bioarrays Files*

Description

Read the header of Codelink files and obtain useful information.

Usage

```
readHeader(file, dec=FALSE)
```

Arguments

file	File to read.
dec	logical; If TRUE determine decimal point.

Value

A list with header and other useful information.

Author(s)

Diego Diez

Examples

```
## Not run:
files <- list.files(pattern = "TXT")
head <- readHeader(files[1])
## End(Not run)
```

readHeaderXLS *Read Header from XLS exported Codelink Bioarrays Files*

Description

Read the header of Codelink files and obtain useful information.

Usage

```
readHeaderXLS(file, dec=FALSE)
```

Arguments

file	File to read.
dec	logical; If TRUE determine decimal point.

Details

This function is not meant to be used by normal users.

Value

A list with header and other usefull information.

Author(s)

Diego Diez

Examples

```
## Not run:
  files <- list.files(pattern = "TXT")
  head <- readHeaderHeader(files[1])
## End(Not run)
```

reportCodelink	<i>Write a report of genes selected in HTML</i>
----------------	---

Description

Takes a list of genes as argument and writes an HTML page containing information about these genes: Unigene, Genbank, Entrez Gene, etc.

Usage

```
reportCodelink(object, chip, filename = NULL, title = "Main title",
               probe.type = FALSE, other = NULL, other.ord = NULL)
```

Arguments

object	an object of class "Codelink" or a list of genes.
chip	the chip description package.
filename	file name used in the report.
title	title used in the report.
probe.type	logical; if TRUE Probe type information is written.
other	list of vectors containing additional values to add to the report.
other.ord	slot name in other to order genes by.

Value

Nothing, only the HTML file generated.

Author(s)

Diego Diez

See Also

[htmlpage](#)

<code>selCV</code>	<i>Select based on CV cutoff.</i>
--------------------	-----------------------------------

Description

Takes a Codelink object and select genes based on CV cutoff.

Usage

```
selCV(object, cutoff)
```

Arguments

<code>object</code>	an object of class "Codelink".
<code>cutoff</code>	cutoff normally calculated with <code>cutCV()</code>

Value

A logical vector.

Author(s)

Diego Diez

<code>SNR</code>	<i>Calculate SNR</i>
------------------	----------------------

Description

Compute SNR inside `read.Codelink`.

Usage

```
SNR(Smean, Bmedian, Bstdev)
```

Arguments

<code>Smean</code>	matrix of Smean intensities.
<code>Bmedian</code>	matrix of background median intensities.
<code>Bstdev</code>	matrix of background standard deviation.

Author(s)

Diego Diez

writeCodelink	<i>Write a Codelink object to file.</i>
---------------	---

Description

Export of the data from a codelink object to a text file.

Usage

```
writeCodelink(object, file, dec = ".", sep = "\t", flag = FALSE, chip)
```

Arguments

object	an object of class "Codelink".
file	filename to write object to.
dec	decimal character to use.
sep	delimiter character to use.
flag	should the Codelink flags be written.
chip	chip package to use, normally guessed.

Details

By default, intensities and SNR are wrote to the file. If set, the flag are also output. The header have "INTENSITY_", "SNR_" and "FLAG_" respectibely appendend to the sample name for those values. The default delimiter is the tab character, but that can be set with the sep argument. The default decimal character is the point.

Author(s)

Diego Diez

Examples

```
## Not run:  
  data(codelink.example)  
  writeCodelink(codelink.example, file = "foo.txt")  
## End(Not run)
```

Index

*Topic **classes**

Codelink-class, 4

*Topic **datasets**

codelink.example, 5

codelink.exprset, 5

*Topic **documentation**

arrayNew, 1

arraySize, 2

as.matrix.Codelink, 2

bkgdCorrect, 3

createWeights, 6

cutCV, 6

decDetect, 7

dim.Codelink, 7

fc2Cond, 8

imageCodelink, 9

logCodelink, 10

mergeArray, 10

na2false, 11

Nondocumented-objects, 12

normalize, 13

normalize.loess, 12

plotCorrelation, 14

plotCV, 15

plotDensities, 15

plotMA, 16

printHead, 17

readCodelink, 18

readHeader, 19

readHeaderXLS, 19

reportCodelink, 20

selCV, 21

SNR, 21

writeCodelink, 22

*Topic **internal**

Nondocumented-objects, 12

*Topic **utilities**

arrayNew, 1

arraySize, 2

as.matrix.Codelink, 2

bkgdCorrect, 3

createWeights, 6

cutCV, 6

decDetect, 7

dim.Codelink, 7

fc2Cond, 8

imageCodelink, 9

logCodelink, 10

mergeArray, 10

na2false, 11

normalize, 13

normalize.loess, 12

plotCorrelation, 14

plotCV, 15

plotDensities, 15

plotMA, 16

printHead, 17

readCodelink, 18

readHeader, 19

readHeaderXLS, 19

reportCodelink, 20

selCV, 21

SNR, 21

writeCodelink, 22

[, Codelink-method

(Codelink-class), 4

annotation, Codelink-method

(Nondocumented-objects), 12

arrayNew, 1

arraySize, 2

as.matrix, 2

as.matrix.Codelink, 2

bkgdCorrect, 3

chipDevice

(Nondocumented-objects), 12

chipDevice, CodelinkSet-method

(Nondocumented-objects), 12

class::CodelinkSet

(Nondocumented-objects), 12

codCorrect

(Nondocumented-objects), 12

codCorrect, CodelinkSet-method

(Nondocumented-objects), 12

Codelink, 5

- Codelink-class, 4
- `codelink.example`, 5
- `codelink.exprset`, 5
- `Codelink2eSet`
(*Nondocumented-objects*), 12
- `CodelinkSet`
(*Nondocumented-objects*), 12
- `CodelinkSet`-class
(*Nondocumented-objects*), 12
- `codNormalize`
(*Nondocumented-objects*), 12
- `codNormalize`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `codPlot` (*Nondocumented-objects*),
12
- `codPlot`, `Codelink`-method
(*Nondocumented-objects*), 12
- `codPlot`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `codPlot`, `MArrayLM`-method
(*Nondocumented-objects*), 12
- `codPreprocess`
(*Nondocumented-objects*), 12
- `codPreprocess`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `codUniverse`
(*Nondocumented-objects*), 12
- `codUniverse`, `Codelink`-method
(*Nondocumented-objects*), 12
- `codUniverse`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `createWeights`, 6
- `cutCV`, 6
- `decDetect`, 7
- `dim.Codelink`, 7
- `eSet2Codelink`
(*Nondocumented-objects*), 12
- `exprs`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `fc2Cond`, 8
- `filterSNR`
(*Nondocumented-objects*), 12
- `filterSNR`, `Codelink`-method
(*Nondocumented-objects*), 12
- `filterSNR`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `filterSNR`, `matrix`-method
(*Nondocumented-objects*), 12
- `getBkg` (*Nondocumented-objects*), 12
- `getBkg`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `getInfo` (*Nondocumented-objects*),
12
- `getInfo`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `getInt` (*Nondocumented-objects*), 12
- `getInt`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `getSNR` (*Nondocumented-objects*), 12
- `getSNR`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `getType` (*Nondocumented-objects*),
12
- `guessAnnotation`
(*Nondocumented-objects*), 12
- `htmlpage`, 20
- `imageCodelink`, 9
- `loess`, 12
- `log2`, 10
- `logCodelink`, 10
- `meanSNR` (*Nondocumented-objects*),
12
- `meanSNR`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `mergeArray`, 10
- `na2false`, 11
- `Nondocumented-objects`, 12
- `normalize`, 13
- `normalize.loess`, 12
- `normalize.median`
(*Nondocumented-objects*), 12
- `plot`, 14
- `plotCorrelation`, 14
- `plotCV`, 15
- `plotDensities`, 15
- `plotDensities`, `Codelink`-method
(*plotDensities*), 15
- `plotMA`, 16
- `printHead`, 17
- `probeNames`
(*Nondocumented-objects*), 12
- `probeNames`, `CodelinkSet`-method
(*Nondocumented-objects*), 12
- `probeTypes`
(*Nondocumented-objects*), 12
- `probeTypes`, `CodelinkSet`-method
(*Nondocumented-objects*), 12

read.table, [18](#)
readCodelink, [18](#)
readCodelinkFiles
 (*Nondocumented-objects*), [12](#)
readCodelinkSet
 (*Nondocumented-objects*), [12](#)
readHeader, [19](#)
readHeaderXLS, [19](#)
reportCodelink, [20](#)

selCV, [21](#)
show, Codelink-method
 (*Codelink-class*), [4](#)
SNR, [21](#)
summary, Codelink-method
 (*Nondocumented-objects*), [12](#)
summaryFlag
 (*Nondocumented-objects*), [12](#)
summaryFlag, Codelink-method
 (*Nondocumented-objects*), [12](#)
summaryFlag, matrix-method
 (*Nondocumented-objects*), [12](#)

writeCluster
 (*Nondocumented-objects*), [12](#)
writeCluster, Codelink-method
 (*Nondocumented-objects*), [12](#)
writeCluster, CodelinkSet-method
 (*Nondocumented-objects*), [12](#)
writeCluster, ExpressionSet-method
 (*Nondocumented-objects*), [12](#)
writeCluster, MArrayLM-method
 (*Nondocumented-objects*), [12](#)
writeCodelink, [22](#)
writeCodelink, CodelinkSet-method
 (*writeCodelink*), [22](#)